





VARILITE®

A non-adjustable skin protection cushion should perform as well after extended use as it does when new. The Reflex shows consistent pressure distribution after rigorous life cycle testing.

FSA pressure map of a new Reflex™ cushion



FSA pressure map of a life-cycled* Reflex™ cushion



*Cycled compressed 20,000 times @ 175lbs

Subject: Sex: Male Height: 6' Weight: 195 lbs Status: Able-bodied

Pressure mapping performed with Force Sensitive Applications (FSA) seat system equipment on July 7, 2007.



4000 1st Avenue South
Seattle, WA 98134
Phone: (800) 827-4548
Fax: (206) 343-5795
www.varilite.com

Authorized EU Representative:
Cascade Designs, Limited
Dwyer Road, Midleton, Co. Cork
Republic of Ireland

