

The Zoid™ Cushion



The Zoid™ Cushion

- Ideal for active users with moderate skin protection needs and mild positioning
- Low profile (2.9 inches)
- Tapered
 - Extra firm foam for lateral stability
 - No medial thigh separator
- 4 way stretch fabric
- Mesh or Incontinence cover
- HCPCS Code E2622 Skin Protection Wheelchair Seat Cushion, adjustable




The Zoid™ Cushion

- Named for its **trapezoid** shape
- 1" (2.5 cm) taper front to back for tighter fit in sports chairs
- 2.9" (7.25 cm) height for a low profile




Zoid™ Foam Stiffness

- Postural support is created by the 3 types of foam:
 - Less firm foam under ischium
 - Medium foam under thighs
 - Firm foam for lateral thigh support



Thigh Area Foam

- Firmer than ischial area foam
- Creates pre-ischial ridge and pelvic bucketing
- Maintains pelvic position by restricting forward and lateral pelvic migration
- Provides support during functional transfers



Lateral Support Foam

- Firmest foam area
- Beveled shape creates bucketing effect of lower extremities to assist sitting balance
- Lateral thigh support minimizes lower extremity abduction
- Lateral pelvic support provides stability
- Postural support created through immersion



The Zoid™ Cover

- Available with mesh or incontinence cover
- Sporty grey fabric
- Cover is designed to protect the thighs of users that sit with legs under the chair.



Clients to Consider:

- Active users with and without sports specific chairs
- Pediatric clients
- Users with knee flexion contractures or who prefer a sitting posture with LEs tucked underneath

