



# Evolution Back™

The back system of choice for wheelchair users with moderate to maximum support needs.



BACK SYSTEMS

## Evolution Back™

Pin bars for VARILITE Evolution Back systems are available in regular or extended lengths. Extended pin bars allow for even greater seat depth and angle adjustability. Extended pin bars are standard with the Evolution Back Deep and optional on Evolution Back regular and tall sizes.



Therapists recommend the Evolution Back because it is easy to mount and remove with one hand, even for the user without fine motor skills. The system is lightweight to conserve energy expenditure for the client who manually propels. Mounting hardware, which comes pre-installed, provides a wide range of lateral, seat-to-back angle, and seat-depth adjustment.




Now Featuring Metal Hardware!




### Upper Pin Bars

#### Evolution Back


##### Upper Pin Bars




Lower Pin Bars




NEUTRAL




POSTERIOR




POSTERIOR



ANTERIOR



NEUTRAL



POSTERIOR

Front of Chair →

Users and caregivers value the Evolution Back's comfort and effective support. Sleek and strong. Lightweight and user friendly.



4000 1st Avenue South  
Seattle, WA 98134  
Phone: (800) 827-4548  
Fax: (206) 343-5795  
Toll Free Fax: (800) 363-7739  
www.varilite.com

Authorized EU Representative:  
Cascade Designs, Limited  
Dwyer Road, Midleton, Co. Cork  
Republic of Ireland

© 2005 Cascade Designs, Inc.



